Quick Guide Dismissing Notifications





Open your email notification and click **"Dismiss**" which will lead you to login to the platform.

38.1 °C	
Thu, 2 Apr 2020 09:08:20 -0400	TB Incubator - Lakewood
	Lab
	Lakewood
	400 New Hampshire Ave, Lakewood, NJ 08701,
	US
	Time Zone America/New_York
	\cap
	Lakewood Incubator is Outside of Safe
	Parameters
	DISMISS





Once clicked, this will take you to the "**Notification List Screen**" where you will find all recent notifications.

On the Home Screen select the

"New Notifications Today" tile.



On the left side, if you select ">" a panel will expand to show you the chart associated to the equipment so you can diagnose the solution.

